

# ORDER FOR CUSTOM-MADE PRODUCTS

## GenuTrain<sup>®</sup>



**F** +44 (0) 121 446-5454  
**E** info@bauerfeind.co.uk

Identifier \_\_\_\_\_  
(e.g. patient name,  
abbreviation) \_\_\_\_\_

**Quotation**  
 **Order**

	Circumference	Length
	f <sup>1</sup>	e - f <sup>1</sup>
	f	e - f
	e	
	d	e - d
	c	e - c

**Knee**  
 Right  Left

**Color**  
 Titan  Nature

**Extras**  
 Siliconized edge\* (add 5 cm to the length)

Tightening strap\*

Velcro fastening\*  
 Top\*\*  Bottom

**Patient details\*\*\***  
 Male  Female  
 Obese  Pronounced muscles

Age \_\_\_\_\_ Height \_\_\_\_\_

\* Subject to additional charges

\*\* Not possible in combination with a siliconized edge; not circular

\*\*\* Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

### All measurements must be entered to ensure proper treatment.

f<sup>1</sup> Thigh: Measure the circumference at least 14 cm above the middle of the knee (e). If the circumference at point f<sup>1</sup> is more than 47 cm, then the circumference must be measured at a higher point, where e - f<sup>1</sup> is at least 17 cm.

f Base of thigh

e Middle of knee

d Narrowest point below the knee

c Calf circumference: Measure the circumference at least 12 cm below the middle of the knee (e).

### Notes:

- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg. The distance between e and f<sup>1</sup> should not be shorter than the distance between e and c.
- If the circumference of f<sup>1</sup> exceeds 60 cm, we recommend attaching a silicone dotted band.

**The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned.**



Stamp

Customer no. \_\_\_\_\_

Signature \_\_\_\_\_

Date \_\_\_\_\_