

ORDERING CUSTOM-MADE PRODUCTS

EpiTrain®




BAUERFEIND.COM

Fax UK ????????

Patient name _____

Quotation only



Arm

right left

Color

titanium natural

Extras

Top band* (+3 cm to the length)

Velcro fastening*

top** bottom

Patient details

male female

obese muscular

Age _____ Height _____

* subject to additional charge

** not possible in combination with siliconized edge; not circular

The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned.

All measurements must be entered to ensure the right treatment is provided.

f¹ Upper arm
(at least 10 cm above e)

f Upper arm sleeve

e Crease of the elbow

d Middle of the forearm muscle

c Muscle attachment
(at least 12 cm below e)

Caution:

- Take measurements with a slight flexion (30°).
- The length is measured on the inside of the arm. Starting point is the crease of the elbow (e).



Stamp

Customer No. _____

Signature _____

Date _____

