

# ORDERING CUSTOM-MADE PRODUCTS

## GenuTrain®



BAUERFEIND.COM

Fax UK ????????

Patient name \_\_\_\_\_  
\_\_\_\_\_

**Quotation only**

	Circumference	Length
	f <sup>1</sup>	e - f <sup>1</sup>
	f	e - f
	e	
	d	e - d
	c	e - c

**Knee**

right  left

**Color**

titanium  natural

**Extras**

Top band\* (+5 cm to the length)

Tightening strap\*

Velcro fastening\*

top\*\*  bottom

**Patient details**

male  female

obese  muscular

Age \_\_\_\_\_ Height \_\_\_\_\_

\* subject to additional charge

\*\* not possible in combination with siliconized edge; not circular

**All measurements must be entered to ensure the right treatment is provided.**

f<sup>1</sup> Thigh  
Measure the circumference at least 14 cm above the middle of the knee (e). If the circumference at point f<sup>1</sup> is more than 47 cm, then the circumference must be measured at ¾ of the length up from the middle of the knee (e).

f Thigh sleeve

e Middle of the knee

d Thinnest point below the knee

c Largest calf circumference  
(at least 12 cm below the middle of the knee)

**Caution:**

- Length measurements are taken with a slight flexion (30°).
- The length is measured on the outside of the leg. The distance between e and f<sup>1</sup> should not be less than the distance between e and c.



Stamp

Customer No. \_\_\_\_\_

Signature \_\_\_\_\_

Date \_\_\_\_\_

