


EFFECTIVE ADJUSTABILITY HOW DOES SecuTec® OA COMPARE?

PROFESSIONAL THERAPY FOR EVERY LEVEL OF TREATMENT

THE PROBLEM WITH CONVENTIONAL ORTHOSES



Front view of right leg

01. DOES NOT ALLOW A PRECISE FIT TO THE LEG

02. THE SPACE BETWEEN ORTHOTIC JOINTS IS NARROWED putting pressure on the painful area

SecuTec® OA SOLVES THE PROBLEM




Front view of right leg

01. FITS TO THE LEG SHAPE

02. TARGETED RELIEF
Orthotic joints retain the space between them to prevent additional pressure on the painful area.



LEVEL I	LEVEL I-II	LEVEL II-(III)	LEVEL III-(IV)	LEVEL IV
GenuTrain A3 mild osteoarthritis of the knee without leg misalignment	GenuTrain S mild osteoarthritis of the knee with joint instability	SofTec OA moderate osteoarthritis of the knee with varus misalignment	SecuTec OA moderate to severe osteoarthritis of the knee with varus or valgus misalignment	MOS-Genu severe osteoarthritis of the knee and after surgery



STANDARD SIZES

Size	1	2	3	4	5	6	7
Circumference A (cm)	38 – 42	41 – 45	44 – 49	48 – 53	51 – 57	54 – 60	60 – 72
Circumference B (cm)	28 – 32	31 – 35	34 – 38	36 – 40	39 – 43	42 – 46	46 – 51

FOR A PERFECT FIT
You can mix and match different frame sizes for upper and lower leg

Size (prox. frame)	1	2	3	4	5	6	7
Size (distal frame)	1,2,3	1,2,3,4	1,2,3,4,5	2,3,4,5,6	3,4,5,6,7	4,5,6,7	5,6,7

ITEM NUMBERS	RIGHT	LEFT
SecuTec OA	121469100506°	121469000507°
Assembly kit	121469200506^>	121469200507^>

°=size, ^... size (proximal), >... size (distal)

<p>BAUERFEIND AG Triebeser Straße 16 07937 Zeulenroda-Triebes Germany P +49 (0) 36628 66-10 00 F +49 (0) 36628 66-19 99 info@bauerfeind.com</p>	<p>UNITED KINGDOM Bauerfeind UK Regents Place 85 Tottenham Court Road London W1T 4TQ P +44 (0) 121 446-5353 F +44 (0) 121 446-5454 info@bauerfeind.co.uk</p>	<p>USA Bauerfeind USA, Inc. 1230 Peachtree St NE Suite 1000 Atlanta, GA 30309 P 1 800 423-3405 F +1 404 201-7800 info@bauerfeindusa.com</p>
--	---	---



I'VE GOT SOMETHING AGAINST OSTEOARTHRITIS

SecuTec® OA
the orthosis for varus and valgus osteoarthritis of the knee

SecuTec® OA

THE ORTHOSIS FOR VARUS AND VALGUS OSTEOARTHRITIS OF THE KNEE

Reduces pain by gently applying pressure to the lower leg and maintaining a constant distance between the orthotic joints.



INDICATIONS

- medial and lateral osteoarthritis of the knee
- postoperative stabilization and relief

TOOL-FREE:

- Extension 0°, 10°, 20°, 30°, 45°
- Flexion 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Fixation 0°, 10°, 20°, 30°, 45° (optional – not included with delivery)



MORE COMFORT FOR YOUR PATIENTS

OPTIMUM PRESSURE DISTRIBUTION

- 01**
- Gel pad can be adapted to your body
 - Additional wedges to ensure a flat fit
- 02**
- Side pad with recess for the head of the fibula



HIGH WEARING COMFORT

- Anti-shock technology in the joints provides a more fluid walking movement
- Microfiber cushioning
- Lies flat against the leg




PERFECT FIT

- Adaptable calf shell tilts and adjusts automatically
- Mix and match frame sizes for upper and lower leg
- Adjustable frame for lower leg



SUCCESSFUL THERAPY



REDUCED PAIN

perceived pain reduced by 36% after 2 weeks of therapy with the SecuTec® OA orthosis¹.

WALK LONGER, BE ACTIVE LONGER



The result of our wearing test: all test subjects who could not walk more than 500 m without an orthosis were able to cover a significantly longer distance after two weeks of therapy with an orthosis.

"IN SPITE OF MY OSTEOARTHRITIS, I'M ABLE TO ENJOY EVERYDAY LIFE AND MY FREE TIME AGAIN."

Learn more about Simon's story at: blog.bauerfeind.com



¹ VAS score reduced by 1.6 points. Unpublished observational study by Bauerfeind AG