

ORDER FOR CUSTOM-MADE PRODUCTS

EpiTrain[®]




F +44 (0) 121 446-5454
E info@bauerfeind.co.uk

Identifier _____
 (e.g. patient name,
 abbreviation) _____

Quotation

Order



Arm

Right Left

Color

Titan Nature

Extras

Siliconized edge* (add 3 cm to the length)

Velcro fastening*

Top** Bottom

Patient details***

Male Female

Obese Pronounced muscles

Age _____ Height _____

* Subject to additional charges

** Not possible in combination with a siliconized edge; not circular

*** Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

All measurements must be entered to ensure proper treatment.

- f^1 Upper arm: Measure the circumference at least 10 cm above the crease of the elbow (e).
- f Base of the upper arm
- e Crease of the elbow
- d Middle of the forearm muscle
- c Base of muscle: Measure the circumference at least 12 cm below the crease of the elbow (e).

Notes:

- Take measurements with a slight flexion (30°).
- The starting point to determine the length measurements is the crease of the elbow (e).
- The custom-made solution does not have a strap.

The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned.



Stamp

Customer no. _____

Signature _____

Date _____