

# ORDER FOR CUSTOM-MADE PRODUCTS

## GenuTrain® P3



Identifier \_\_\_\_\_  
(e.g. patient name)

**F** +44 (0) 121 446-5454  
**E** info@bauerfeind.co.uk

**Quotation**       **Order**

**Knee**

Right       Left

**Color**

Titan       Nature

**Extras**

Siliconized edge\* (add 5 cm to the length)

**Position of the corrective guide**

Standard (lateral)       Alternative option (medial)

**Patient details\*\***

Male       Female  
 Obese       Pronounced muscles

Age: \_\_\_\_\_ Height: \_\_\_\_\_ cm

\* Subject to additional charges – please consider the additional length when measuring patients

\*\* Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

The purchaser is responsible for ensuring that all size information is correct.  
Custom-made solutions cannot be returned.

**All measurements must be entered to ensure proper treatment.**

- f<sup>1</sup> Thigh / top edge of the support
- f Base of thigh (approx. midpoint between e and f<sup>1</sup>)
- e Middle of knee
- d Narrowest point below the knee, approx. 5–6 cm
- c (Largest) calf circumference
- b<sup>1</sup> Bottom edge of the support



Stamp

Customer no. \_\_\_\_\_

Signature \_\_\_\_\_

Date \_\_\_\_\_

**short** (max. body height approx. 170 cm)

	Circumferences	Lengths
	f <sup>1</sup>	e – f <sup>1</sup> = 18 cm
	f	e – f
	e	
	d	e – d
	c	e – c
	b <sup>1</sup>	e – b <sup>1</sup> = 15 cm

**long** (min. body height approx. 170 cm)

	Circumferences	Lengths
	f <sup>1</sup>	e – f <sup>1</sup> = 22 cm
	f	e – f
	e	
	d	e – d
	c	e – c
	b <sup>1</sup>	e – b <sup>1</sup> = 17 cm

**Notes:**

- Take measurements while standing, with a slight flexion (30°).
- The lengths are measured on the outside of the leg.

**Please note:**

- **Min. circumf. f<sup>1</sup> = 39 cm**
- **Min. circumf. f = 37 cm**
- **Min. circumf. e = 30 cm**
- **Min. circumf. b<sup>1</sup> = 28 cm**