

# ORDER FOR CUSTOM-MADE PRODUCTS

## GenuTrain<sup>®</sup> S



Identifier \_\_\_\_\_  
(e.g. patient name,  
abbreviation) \_\_\_\_\_

**F** +44 (0) 121 446-5454  
**E** info@bauerfeind.co.uk

**Quotation**

**Order**

	Circumference	Length
	<b>f<sup>1</sup></b>	<b>short: e - f<sup>1</sup> = 18 cm</b> <b>long: e - f<sup>1</sup> = 21 cm</b>
	<b>f</b>	<b>e - f</b>
	<b>e</b>	
	<b>d</b>	<b>e - d</b>
	<b>c</b>	<b>e - c</b>
	<b>b<sup>1</sup></b>	<b>e - b<sup>1</sup> = 18 cm</b>

### Knee

Left  Right

short  
(Length of support, thigh [e - f<sup>1</sup>] = 18 cm)

long  
(Length of support, thigh [e - f<sup>1</sup>] = 21 cm)

### Extras

Siliconized edge\* (add 5 cm to the length)

### Patient details\*\*

Male  Female

Obese  Pronounced muscles

Age \_\_\_\_\_ Height \_\_\_\_\_

\* Subject to additional charges – please consider the additional length when measuring patients

\*\* Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

### All measurements must be entered to ensure proper treatment.

**f<sup>1</sup>** Thigh  
short version: Measure the circumference 18 cm above the middle of the knee (e).  
long version: Measure the circumference 21 cm above the middle of the knee (e).

**f** Base of thigh

**e** Middle of knee

**d** Narrowest point below the knee

**c** Calf circumference: Measure the circumference at least 12 cm below the middle of the knee (e).

**b<sup>1</sup>** Calf circumference 18 cm below the knee

### Notes:

- Custom-made solutions are only available in the color titan.
- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg.
- The selection of the short / long version depends on the individual measurements of the patient. The long version is recommended for patients from a height of 1.80 m.
- The following minimum circumference measurements are required for supplying the product:  
f (14 cm above the knee): at least 38 cm  
e (middle of knee): at least 32.5 cm  
c (12 cm below the knee): at least 28 cm



Stamp

Customer no. \_\_\_\_\_

Signature \_\_\_\_\_

Date \_\_\_\_\_