

ORDER FOR CUSTOM-MADE PRODUCTS

GenuTrain® S Pro



Identifier _____
(e.g. patient name,
abbreviation) _____

F +44 (0) 121 446-5454
E info@bauerfeind.co.uk

Quotation

Order

	Circumference	Length
	f¹	short: e - f¹ = 18 cm long: e - f¹ = 21 cm
	f	e - f
	e	
	d	e - d
	c	e - c
	b¹	e - b¹ = 18 cm

Knee

Left Right

short
(Length of support, thigh [e - f¹] = 18 cm)

long
(Length of support, thigh [e - f¹] = 21 cm)

Extras

Siliconized edge* (add 5 cm to the length)

Patient details**

Male Female

Obese Pronounced muscles

Age _____ Height _____

* Subject to additional charges – please consider the additional length when measuring patients

** Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

All measurements must be entered to ensure proper treatment.

f¹ Thigh
short version: Measure the circumference 18 cm above the middle of the knee (e).
long version: Measure the circumference 21 cm above the middle of the knee (e).

f Base of thigh

e Middle of knee

d Narrowest point below the knee

c Calf circumference: Measure the circumference at least 12 cm below the middle of the knee (e).

b¹ Calf circumference 18 cm below the knee

Notes:

- Custom-made solutions are only available in the color titan.
- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg.
- The selection of the short / long version depends on the individual measurements of the patient. The long version is recommended for patients from a height of 1.80 m.
- The following minimum circumference measurements are required for supplying the product:
f (14 cm above the knee): at least 38 cm
e (middle of knee): at least 32.5 cm
c (12 cm below the knee): at least 28 cm



Stamp

Customer no. _____

Signature _____

Date _____