

ORDER FOR CUSTOM-MADE PRODUCTS

MalleoTrain®



F +44 (0) 121 446-5454
E info@bauerfeind.co.uk

Identifier _____
 (e.g. patient name,
 abbreviation) _____

Quotation

Order

Foot

Right Left

Color

Titan Nature

Extras

Velcro fastening*

Measuring points for pads**

① _____

② _____

Patient details***

Male Female
 Obese Pronounced muscles

Age _____ Height _____

* Subject to additional charges

** Special position of the pads (position of ankle bones differ substantially from the standard measurements)

*** Optional; please specify special anatomical features to ensure an accurate fit of the finished product.

The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned.




Stamp

Customer no. _____

Signature _____

Date _____

Circumference **Length**



All measurements must be entered to ensure proper treatment.

- b¹ Base of calf: Measure the circumference at least 16 cm above the heel (A).
- b Smallest circumference above the ankle bone
- h Instep (circular measurement of the instep over the heel)
- a Circumference behind the MTP joint of the little toe
- A Heel
- A-a Foot length without toe

Notes:

- Only specify the pad measurement points if special positioning is required.
- Distance ① from the medial malleolus to the lateral malleolus (measured horizontally across the Achilles tendon)
- Distance ② from the medial malleolus to the lateral malleolus (measured over the sole of the foot)