

PERFECT FIT SIZING SYSTEM



VenoTrain® Round-knit products

BAUERFEIND.CO.UK

Circumference	XS*		S		M		L		XL	
	Normal	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus	
cG top band, wide		51-57	59-66	56-63	62-69	60-67	66-74	63-71	71-80	
cG top band, standard	44-50	45-51	52-59	49-56	55-62	53-60	58-66	55-63	62-71	
cG edging, AT, top band**		45-57	52-66	49-63	55-69	53-67	58-74	55-71	62-80	
cF	39-50	40-52	43-55	43-55	46-60	45-59	50-64	49-63	52-68	
cE	32-38	33-40	36-44	35-43	39-49	38-48	42-53	42-53	47-58	
cD	28-31	29-35	32-38	32-38	35-42	34-42	37-45	37-45	41-49	
cC	29-33	31-37	34-40	34-40	37-44	37-45	39-47	40-48	43-51	
cB1	23-26	25-29	27-31	27-32	29-34	29-35	31-37	32-38	34-40	
cB	18-20	20-22		22-25		25-28		28-31		
cY	27-31	29-33		31-35		34-40		37-44		
cA	17-22	19-24		22-27		25-31		28-35		

Length	short	long	Type
lG	61-70	70-80	AG, AT
lD	35-40	40-46	AD

* only for VenoTrain soft S – limited product range
 ** Top band lace sensitive for VenoTrain micro Design Edition

VenoTrain micro | VenoTrain micro balance | VenoTrain look | VenoTrain business | VenoTrain soft | VenoTrain soft S

The Perfect-Fit sizing system developed by Bauerfeind on the basis of digital leg measurements offers unique flexibility in ensuring that medical compression stockings from the standard product range fit perfectly on your legs.



All measurements in the sizing system are given in cm: l stands for length · c stands for circumference · Y stands for instep

*For information on the custom-made solutions available, please refer to the current price list.