

V6+ SIZING SYSTEM AND PRODUCT RANGE

VenoTrain® impuls

Circumference	I		II		III		IV		V		VI	
	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus
cG	45-53	51-59	48-56	54-62	51-59	57-66	54-63	60-69	57-66	63-73	60-70	66-76
cF	39-45	43-50	42-49	46-53	45-52	49-57	48-56	52-60	51-59	55-64	54-63	58-67
cE	32-37	35-40	34-39	37-42	36-41	39-44	38-44	41-47	40-46	43-49	42-48	45-51
cD	28-32	30-34	30-34	32-37	32-37	34-39	34-39	36-41	36-41	38-43	38-43	40-46
cC	29-33	31-35	31-35	33-38	33-38	35-40	35-40	37-42	37-42	39-44	39-44	41-47
cB1	23-27	24-28	25-29	26-30	27-31	28-32	29-34	30-35	31-36	32-37	33-38	34-39
cB	19-20.5		20.5-22		22-24		24-26		26-28.5		28.5-31	
cY	25-29		26-31		28-33		30-35		32-38		35-41	
cA	19-22		20-23		22-25		24-27		26-30		28-33	

Length	short	long	Type
ℓG	62-70	71-80	AG
ℓD	35-39	40-44	AD

Effect	Properties	Compression classes	Design*	Toe	Edging/top band	Color
Increased compression For advanced vein problems and mild lymphatic disorders	<ul style="list-style-type: none"> • Quicker pressure increase for high working pressure • Comfortably soft and durable • Highly differentiated standard sizing system 	Ccl 2 3	AD AG	Open	Top band Silicone dots Garter belt (AG)	CARAMEL BLACK

Express manufacture

All measurements in the sizing system are given in cm: ℓ stands for length · c stands for circumference · Y stands for instep

*For information on the custom-made solutions available, please refer to the current price list.