


# SIZING SYSTEM AND PRODUCT RANGE

## VenoTrain® delight

BAUERFEIND.CO.UK



Circumference	I	II	III	IV	V	VI
<b>cD</b>	28–34	30–37	32–39	34–41	36–43	38–46
<b>cC</b>	29–35	31–38	33–40	35–42	37–44	39–47
<b>cB1</b>	23–28	25–30	27–32	29–35	31–37	33–39
<b>cB</b>	19–20,5	20,5–22	22–24	24–26	26–28,5	28,5–31
<b>cY</b>	25–29	26–31	28–33	30–35	32–38	35–41
<b>cA</b>	19–22	20–23	22–25	24–27	26–30	28–33

Length	short	long	Type
ℓD	35 – 39	40 – 44	AD

Effect	Properties	Compression classes	Design*	Toe	Edging/top band	Color
<b>Increased compression</b> Compression stocking for intensive therapy	<ul style="list-style-type: none"> <li>Breathable double-stretch knit</li> <li>Latex-free and gentle on the skin</li> <li>Extra-flat seams</li> </ul>	Ccl 2   3	AD	Open Closed	Top band Silicone dots	<div style="display: flex; gap: 5px;"> <div style="background-color: #ccc; border-radius: 10px; padding: 2px 10px; display: inline-block;">CAMEL</div> <div style="background-color: #ccc; border-radius: 10px; padding: 2px 10px; display: inline-block;">BLACK</div> </div>