


X5+ SIZING SYSTEM AND PRODUCT RANGE

VenoTrain® ulcertec

BAUERFEIND.CO.UK



Circumference	XS		S		M		L		XL	
	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus
cC	24-29	28-34	27-32	31-38	29-36	35-42	33-40	38-47	36-44	43-52
cB	20-21.5		21.5-23		23-25		25-28		28-32	
Length	short		long		Type					
ℓD	35-40		40-45		AD					

Effect	Properties	Compression classes	Design	Toe	Edging	Color
Venous ulcers For healed and active venous ulcers	Compression stocking set: Outer stocking (patented rhomboid knit) Liner (skin-friendly material) Medical effectiveness proven by clinical study	VenoTrain® ulcertec 39 – moderate Compression strength 34-46 mmHg, liner with gentle compression • Ulcers with a good tendency to heal • Light liner for patients who only need little compression at night	AD	Open outer stocking	Top band	NATURAL
				Closed liner	Top band	WHITE
		VenoTrain® ulcertec 46 – strong Compression strength 34-46 mmHg, liner with moderate compression • Ulcers with a poor tendency to heal • Liner with moderate compression for those able to tolerate greater compression – including at night				